



INVACARE® TOP END®
EXCELERATOR™ SERIES
HANDCYCLE



If you want a great way to exercise, cross-train or just have fun, one of the **Invacare® Top End® Excelerator™ Series Handcycles** is just what you need. This upright and extremely maneuverable three-wheel handcycle is capable of speeds up to 15 mph. The smooth seven-speed Shimano® hub with its reversing drum brake delivers hands-on braking and allows the rider to tackle hills or level terrain with equal ease. The oval tube carbon steel frame construction with adjustable sliding seat makes transfers in and out easier too.

The Excelerator Handcycle is stocked in three seat widths: 16", 18" or 20" wide, plus you can order the popular mountain drive option which adds 7 extra low gears for climbing hills. No more long lead times so...Get Cranking!

Features

Carbon steel frame and fork
Adjustable seat: forwards and backwards
Adjustable back height
Fold over back upholstery
Adjustable crank height
Chain tension idler
Full chain guard, reflectors and safety flag
Adjustable height/angle footrest system
Shimano® 7 speed hub w/reversing drum brake

Parking brake
175 mm crankset w/one chain ring
Ovalized ergonomic hand pedals
Performance wheels w/cruiser tires
Quick release axles
3° camber
Cushion and hook and loop seat restraint
Weight capacity 350 lb.

Model #: XCL



1. Twist Shifter for 7 speed hub
2. Quad Twist shifter option
3. Parking brake
4. Mountain drive transmission option
5. Frictionless ovalized handpedal option



Options

Mountain drive transmission kit adds 7 low gears
Quad twist shifter helps those with hand impairments
Tri-pin quad handpedals
Crank width adapter set
Quad gloves
Water bottle and cage
Tow bar for wheelchair
Cordless computer

Front safety light
Rear flashing light
Helmet
Seat restraint: Auto style fastener
Mirror
Crutch holder and strap
Bike rack for trailer hitch
Top End t-shirt

Features

Carbon steel frame and fork

Adjustable seat: forwards and backwards

Adjustable back height

Fold over back upholstery

Adjustable crank height

Chain tension idler

Full chain guard, reflectors and safety flag

Adjustable height/angle footrest system

Shimano® 7 speed hub with reversing drum brake

Parking brake

155 mm crankset with one chain ring

Foam covered horizontal handpedals

20" Cruiser wheels and tires

Quick release axles

9° camber

Cushion and hook and loop seat restraint

Weight capacity 250 lb.

Model #: LXC



1. Parking brake

2. Ergonomic handpedals

3. Hook and loop safety seat restraint



Options

Quad twist shifter option

Quad gloves

Mirror

Water bottle and cage

Helmet

Front safety light

Rear flashing light

Crutch holder and strap

Bike rack for trailer hitch

Seat restraint: Auto style fastener

Features



1. V Crankset

The Invacare® Top End® V crankset option accommodates consumers who might have “clearance” issues. Not compatible with mountain drive option



2. Mountain Drive Transmission

Available on the Excelerator handcycle
Adds 7 low gears (7 low + 7 standard = 14 total)
Easy shifting, especially for quadriplegics – just push-in the center button
Lightweight planetary transmission
A must have for hilly terrain



3. Full Chain Guard

A full chain guard on all 7 speed handcycles is a great safety feature

Invacare Corporation
www.invacare.com

USA

One Invacare Way
Elyria, Ohio
44035-4190
(800) 333-6900

Canada

570 Matheson Blvd. E.,
Unit 8
Mississauga, Ontario
L4Z 4G4 Canada
(800) 668-5324

Top End

www.topendwheelchair.com
Sports and Recreation
Products
4501 63rd Circle North
Pinellas Park, Florida
33781-5914
(727) 522-8677
(800) 532-8677
Fax: (727) 522-1007

©2013

Invacare Corporation.
All rights reserved.
Trademarks are identified
by the symbols ™ and ®.
All trademarks are owned
by or licensed to Invacare
Corporation unless
otherwise noted. Shimano
is a registered trademark of
Shimano Inc.
Form No. 12-509
120736

